

The Thermo Scientific Orion 3-Star Plus pH meters are ideal for a variety of pH measurements – water, media, beverage, process water.

## Thermo Scientific Orion 3-Star Plus

pH Meters



### Features and Benefits

- Simultaneous display of pH or mV, and temperature measurements on the backlit LCD
- The Thermo Scientific Orion SMART STABILITY™ and SMART AVERAGING™ functions for automatically optimized accuracy, precision and response time
- Temperature calibration capability
- pH calibration of up to 5 points with calibration alarm
- Automatic buffer recognition of USA, NIST and DIN buffers
- Storage of up to 10 individually password protected methods for easy retrieval of operation procedures
- Datalog up to 1000 points with time and date stamp
- Average slope and slope of segment printed with pH measurements
- RS232 port for data downloads and communication
- Benchtop units can control an autosampler and the 096010 stirrer probe (*each sold separately*)
- Benchtop units are splashproof with an IP54-rated housing and include a universal power supply
- Portable units are waterproof with an IP67-rated housing and run for over 2,000 hours on four AA batteries
- 3 year meter warranty

## Product Specifications and Ordering Information


### Thermo Scientific Orion 3-Star Plus pH Meters

<b>pH</b>	Range	- 2.000 to 19.999
	Resolution	0.1, 0.01, 0.001
	Relative Accuracy	± 0.002
<b>mV</b>	Range	± 1999.9
	Resolution	0.1
	Relative Accuracy	± 0.2 or 0.05%, whichever is greater
<b>Temperature</b>	Range	- 5 to 105 °C
	Resolution	0.1 up to 99.9 °C, 1.0 over 99.9 °C
	Relative Accuracy	± 0.1 °C

### Ordering Information

Meters Only		
Cat. No.	Package	
1112000	Includes benchtop meter, universal power and user guide	
1212000	Includes portable meter, batteries and user guide	
Portable Meter Kits		
Cat. No. / Application	Sensor (s)	Accessories and Solutions
1212001 / Fresh Water	<ul style="list-style-type: none"> <li>9107BNMD Low maintenance epoxy pH/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Beaker</li> <li>Electrode storage chamber</li> <li>pH electrode storage solution (5 x 60 mL bottles)</li> <li>pH 4.01, 7.00, 10.01 buffers (10 x 15 mL pouches)</li> <li>Hard field case</li> </ul>
1212002 / TRIS Clog-Free Wastewater	<ul style="list-style-type: none"> <li>5107BNMD No Cal refillable pH/ATC electrode</li> </ul>	<ul style="list-style-type: none"> <li>Beaker</li> <li>pH electrode storage solution (5 x 60 mL bottles)</li> <li>pH 4.01, 7.00, 10.01 buffers (10 x 15 mL pouches)</li> <li>Hard field case</li> </ul>
1212501 / Fresh Water	<ul style="list-style-type: none"> <li>9107BNMD Low maintenance epoxy pH/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Beaker</li> <li>Electrode storage chamber</li> <li>pH electrode storage solution (5 x 60 mL bottles)</li> <li>pH 4.01, 7.00, 10.01 buffers (10 x 15 mL pouches)</li> <li>Soft field case</li> </ul>
1212503 / ORP	<ul style="list-style-type: none"> <li>9179BNMD Low maintenance epoxy ORP/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Meter impact jacket with electrode storage chamber</li> <li>ORP standard (5 x 60 mL bottles)</li> <li>Soft field case</li> </ul>
1212101 / Fresh Water	<ul style="list-style-type: none"> <li>9107BNMD Low maintenance epoxy pH/ATC Triode electrode</li> </ul>	—
1010115 / Fresh Water	<ul style="list-style-type: none"> <li>9107BNMD Low maintenance epoxy pH/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Meter impact jacket with electrode storage chamber</li> <li>Electrode rinse solution (25 x 15 mL pouches)</li> <li>pH electrode storage solution (5 x 60 mL bottles)</li> <li>pH 4.01, 7.00 and 10.01 buffers (10 x 15 mL pouches)</li> <li>Soft field case</li> <li>Method Note - pH of Salt Water</li> </ul>
1010110 / Wastewater	<ul style="list-style-type: none"> <li>5107BNMD No Cal refillable pH/ATC electrode</li> </ul>	<ul style="list-style-type: none"> <li>Electrode storage chamber</li> <li>pH electrode storage solution (5 x 60 mL bottles)</li> <li>pH 4.01, 7.00 and 10.01 buffers (10 x 15 mL pouches)</li> <li>Method Note - pH of Wastewater</li> </ul>

### Ordering Information (cont'd)



#### Benchtop Meter Kits

Cat. No. / Application	Sensor (s)	Accessories and Solutions
1112001 / Fresh Water	<ul style="list-style-type: none"> <li>9157BNMD refillable epoxy pH/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> <li>pH 4.01, 7.00, and 10.01 buffers and pH electrode storage solution (pint bottles)</li> </ul>
1112003 / TRIS Fresh Water	<ul style="list-style-type: none"> <li>8102BNUWP ROSS Ultra refillable glass pH electrode</li> <li>927005MD Epoxy ATC probe</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> <li>pH 4.01, 7.00, and 10.01 buffers and pH electrode storage solution (pint bottles)</li> </ul>
1112004 / Clog-Free Wastewater	<ul style="list-style-type: none"> <li>9172BNWP Sure-Flow refillable glass pH electrode</li> <li>927007MD Stainless steel ATC probe</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> <li>pH 4.01, 7.00, and 10.01 buffers and pH electrode storage solution (pint bottles)</li> </ul>
1112006 / TRIS Ultra Pure Water	<ul style="list-style-type: none"> <li>8102BNUWP ROSS Ultra refillable glass pH electrode</li> <li>927005MD Epoxy ATC probe</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> <li>pH 4.10, 6.79, and 9.15 pure water buffers and pH electrode storage solution (pint bottles)</li> <li>pHISA adjuster solution (5 x 60 mL bottles)</li> </ul>
1112201 / Fresh Water	<ul style="list-style-type: none"> <li>9157BNMD Refillable epoxy pH/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> </ul>
1112202 / Clog-Free TRIS Wastewater	<ul style="list-style-type: none"> <li>8172BNWP ROSS Sure-Flow refillable glass pH electrode</li> <li>927006MD Glass ATC probe</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> </ul>
1112203 / TRIS Fresh Water	<ul style="list-style-type: none"> <li>9156APWP AquaPro low maintenance epoxy pH electrode</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> </ul>
1112204 / Fresh Water	<ul style="list-style-type: none"> <li>9107BNMD Low maintenance epoxy pH/ATC Triode electrode</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> </ul>
1112205 / Fresh Water	<ul style="list-style-type: none"> <li>9156BNWP Refillable epoxy pH electrode</li> <li>927007MD Stainless steel ATC probe</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> </ul>
1112207 / Fresh Water	<ul style="list-style-type: none"> <li>9106BNWP Low maintenance epoxy pH electrode</li> </ul>	<ul style="list-style-type: none"> <li>Swing arm electrode stand</li> </ul>
1112101 / TRIS Fresh Water	<ul style="list-style-type: none"> <li>8102BNUWP ROSS Ultra refillable glass pH electrode</li> <li>927007MD Stainless steel ATC probe</li> </ul>	—
1112102 / Clog-Free TRIS Wastewater	<ul style="list-style-type: none"> <li>8172BNWP ROSS Sure-Flow refillable glass pH electrode</li> <li>927006MD Glass ATC probe</li> </ul>	—
1112103 / TRIS Fresh Water	<ul style="list-style-type: none"> <li>8102BNUWP ROSS Ultra refillable glass pH electrode</li> </ul>	—
1112104 / Fresh Water	<ul style="list-style-type: none"> <li>9157BNMD Refillable epoxy pH/ATC Triode electrode</li> </ul>	—
1112107 / Fresh Water	<ul style="list-style-type: none"> <li>9106BNWP Low maintenance epoxy pH electrode</li> </ul>	—
1112108 / Clog-Free Wastewater	<ul style="list-style-type: none"> <li>9165BNWP Sure-Flow refillable epoxy pH electrode</li> </ul>	—

©2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. ROSS and the CDIL tradenames are trademarks of Thermo Fisher Scientific Inc. US Patent 6,793,787.



S-3STARPH-E 0408 RevA

Environmental Instruments  
Water Analysis Instruments

**North America**  
166 Cummings Center  
Beverly, MA 01915 USA  
Toll Free: 1-800-225-1480  
Tel: 1-978-232-6000  
Dom. Fax: 1-978-232-6015  
Int'l Fax: 978-232-6031

**Europe**  
Denmark House, Angel Drive  
Ely, Cambridgeshire  
England, CB7 4ET  
Tel: 44-1353-666111  
Fax: 44-1353-666001

**Asia Pacific**  
Blk 55, Ayer Rajah Crescent  
#04-16/24, Singapore 139949  
Tel: 65-6778-6876  
Fax: 65-6773-0836

[www.thermo.com/water](http://www.thermo.com/water)

**Thermo**  
SCIENTIFIC